


Instrument Trays

Minitainer Cases

The Liberty standard cases portfolio offers you cost effective solutions for organising and protecting your valuable instruments. The range of cases and trays are made of Radel® R, aluminium or stainless steel. In addition, these cases have FDA 510(k) clearance for use in gravity or pre-vacuum steam and ETO sterilisation and are compatible with CSR wrap and filtered containers.

The Liberty range of cases and trays utilise mats and bars, and will protect your micro-surgical, ENT, orthopedic or other instruments that may otherwise become damaged if left loose in wire baskets or containers.

One of these systems offered by Liberty is the range of Universal Cases. These cases are designed with trays to accommodate your micro to medium sized instruments. The universal trays provide ample space to configure comprehensive procedural sets. These cases are ideally suited for use with neurosurgical, ear/nose/throat, reconstructive, trauma and arthroscopic instruments.



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: LSB1-Instrument-Trays_v4-1

First Prepared: 11/03/15 Eevien Tan
Release Approval: 11/03/15 Lee Carter
Version Control: V4.1 Page 1 | 3

Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice.
Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE
Visit our website at: <http://www.lrinstruments.com.au>

Instrument Trays

Minitainer Cases



LSB1-2700 Minitainer, 2 x 7 Base, Lid and Mat
DIM Internal 181 x 49 x 18.5 MM



LSB1-3700 Minitainer, 3x7 Base, Lid and Mat
DIM Internal 180 x 73 x 18.5 MM



LSB1-2710 Minitainer, 2x7 Base with 2 Silicone Strips and Lid
DIM Internal 181 x 49 x 18.5 MM



LSB1-3710 Minitainer, 3x7 Base with 2 Silicone Strips and Lid for 3 Diamond Knives
DIM Internal 180 x 73 x 18.5 MM



LSB1-2720 Minitainer, 2x7 Base with 2 Silicone Strips and Lid for 2 Diamond Knives
DIM Internal 181 x 49 x 18.5 MM



LSB1-3720 Minitainer, 3x7 Base with 2 Silicone Strips and Lid for 9 Cannulae
DIM Internal 180 x 73 x 18.5 MM



LSB1-3600 Minitainer, 3x6 Base, Lid and Mat
DIM Internal 148 x 64 x 18 MM



LSB1-4800 Minitainer, 4 x 8 Base, Lid and Mat
DIM Internal 190 x 102 x 19 MM



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: LSB1-Instrument-Trays_v4-1

First Prepared: 11/03/15 Eevien Tan
Release Approval: 11/03/15 Lee Carter
Version Control: V4.1 Page 2 | 3

Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice.
Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE
Visit our website at: <http://www.lrinstruments.com.au>


Instrument Trays

Minitainer Cases



LSB1-7100 Minitainer, 7 x 11 Single Level Base, Lid and Shallow Mat

DIM Internal 267 x 159 x 16 MM



LSB1-7120 Minitainer, 7 x 11 Double Level Base, Lid and 2 Mats

DIM Internal 267 x 159 x 35 MM



LSB1-7110 Minitainer, 7 x 11 Single Level Base, Lid and Deep Mat

DIM Internal 267 X 159 X 35 MM

Code	Description
LSB1-2700	POLYVAC MINITAINER 181 X 49 X 18.5MM (2 X 7) BASE, LID AND MAT
LSB1-2710	POLYVAC MINITAINER 181 X 49 X 18.5MM (2 X 7) BASE WITH 2 SILICONE STRIPS AND LID FOR 1 DIAMOND KNIFE
LSB1-2720	POLYVAC MINITAINER 181 X 49 X 18.5MM (2 X 7) BASE WITH 2 SILICONE STRIPS AND LID FOR 2 DIAMOND KNIVES
LSB1-3600	POLYVAC MINITAINER 148 X 64 X 18MM (3 X 6) BASE, LID AND MAT
LSB1-3700	POLYVAC MINITAINER 180 X 73 X 18.5MM (3 X 7) BASE, LID AND MAT
LSB1-3710	POLYVAC MINITAINER 180 X 73 X 18.5MM (3 X 7) BASE WITH 2 SILICONE STRIPS AND LID FOR 3 DIAMOND KNIVES
LSB1-3720	POLYVAC MINITAINER 180 X 73 X 18.5MM (3 X 7) BASE WITH 2 SILICONE STRIPS AND LID FOR 8 CANNULAE
LSB1-4800	POLYVAC MINITAINER 190 X 102 X 19MM (4 X 8) BASE LID AND MAT
LSB1-7100	POLYVAC MINITAINER 267 X 159 X 16MM (7 X 11) SINGLE LEVEL BASE, LID AND SHALLOW MAT
LSB1-7110	POLYVAC MINITAINER 267 X 159 X 35MM (7 X 11) SINGLE LEVEL BASE, LID AND DEEP MAT
LSB1-7120	POLYVAC MINITAINER 267 X 159 X 35MM (7 X 11) DOUBLE LEVEL BASE, LID AND 2 MATS



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: LSB1-Instrument-Trays_v4-1

First Prepared: 11/03/15 Eevien Tan
Release Approval: 11/03/15 Lee Carter
Version Control: V4.1 Page 3 | 3

Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice.
Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE
Visit our website at: <http://www.lrinstruments.com.au>