

# Cardiology Ultrasharp Liberty Stethoscope



## Professional Cardiology Stethoscope

- Ergonomic and unique design
- Two-in-one tube system for improved sound conduction
- Comfortable and flexible binaural
- Silicone ear tips make the auscultation clearer

## Appearance Features

Concise, sleek and smooth ergonomic design for that professional touch.

## Structure Features

Special Zinc Alloy chest piece is beneficial for sound collecting and conducting. Electroplating technology make the surface more smooth and effectively resists tarnishing.

The length of spring is 60mm in it's natural state and comfortable to wear. The angle of the ear tube is set at 75°.

The audio-conductivity is particularly good. Extra strong A-tube made of PVC have separate sound conducting pathways for the left and right ears.

Silicone fitted ear tip, sounds more clearly, is durable and comfortable to wear. The internal benched design can avoid the ear tube penetrating the ear tip.



## Technical Details

Material	Head Piece	Electroplated
	Chest Piece	Zinc Alloy
	Tubing	Silicone, Latex Free
	Ear Tips	PVC
	Diaphragm	Hypo-allergenic insulating polymer
Dimensions	Length	80cm
	Weight	Complete Unit 106gms
AU Distributor	Liberty Health	
Warranty	1 Year Whole Set	
	5 Years Chest Piece	
	5 Years Ear Tube	



Chest Piece



Ear Tips

Product Code	Description
LSCAB	CARDIOLOGY ULTRASHARP STETHOSCOPE BLACK

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134  
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au  
File Name: LSCAB\_Cardiology-Ultrasharp\_v2-0



First Prepared: 21/01/15 Eevien Tan  
Release Approval: 21/01/15 Lee Carter  
Version Control: V2.0 Page 1 | 1