

Classic Turnable Stethoscope

With Pressure Sensitive Diaphragm



Dual head professional stethoscope with superior acoustic sensitivity and pressure sensitive diaphragm. The high performance stainless steel chest piece is at the top of its class. The definitive doctors stethoscope.

New and Innovative Stethoscope

Light non-corrosive titanium design features "Superflex" binaural tubing assembly, colour coordinated plastic retaining, and non-chill rings. Complete with the innovate pressure sensitive diaphragm and advanced ear tip.

Pressure Sensitive Diaphragm

This new design offers a pressure sensitive diagnosis, allowing clearer diagnosis of both high and low frequency sounds without having to change sides, as with standard stethoscopes.

Advanced Soft Eartip

With its patented two-layered design, this model is extra flexible for added comfort.

Binaural

"Superflex" deluxe stainless inner spring binaural assembly with ear tip.

Features

- Tunable diaphragm
- Innovative design
- Optimum acoustic performance
- Finest craft finish
- Multi patents
- Soft ear tips
- Dual Head Chest piece with non-chill ring
- Latex Free
- 5 year warranty
- Free name tag with each stethoscope
- Free spare ear tips and diaphragm
- Professionally packaged in a cardboard box



Chest Piece



Ear Tips



| Product Code | Description |
|--------------|--|
| LSCLTB | CLASSIC TUNABLE STETHOSCOPE BLACK |
| LSCLTNB | CLASSIC TUNABLE STETHOSCOPE NAVY BLUE |
| LSCLTO | CLASSIC TUNABLE STETHOSCOPE ORANGE |
| LSCLTPI-MF | CLASSIC TUNABLE STETHOSCOPE MCGRATH PINK |
| LSCLTPU | CLASSIC TUNABLE STETHOSCOPE PURPLE |
| LSCLTR | CLASSIC TUNABLE STETHOSCOPE RASPBERRY |
| LSCLTRB | CLASSIC TUNABLE STETHOSCOPE ROYAL BLUE |



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: LSCLT_Classic-Turnable_v2-0

First Prepared: 21/01/15 Eevien Tan
 Release Approval: 21/01/15 Lee Carter
 Version Control: V2.0 Page 1 | 2

Classic Turnable Stethoscope

With Pressure Sensitive Diaphragm



Features

- Tunable diaphragm
- Innovative design
- Optimum acoustic performance
- Finest craft finish
- Multi patents
- Soft ear tips
- Dual Head Chest piece with non-chill ring
- Latex Free
- 5 year warranty
- Free name tag with each stethoscope
- Free spare ear tips and diaphragm
- Professionally packaged in a cardboard box

Technical Details

| | | |
|----------------|----------------|------------------------------------|
| Material | Head Piece | Non-corrosive Titanium |
| | Chest Piece | Stainless steel brush finished |
| | Tubing | PVC |
| | Ear Tips | Latex Free |
| | Diaphragm | Hypo-allergenic insulating polymer |
| Dimensions | Length | 77cm |
| Weight | Complete Unit | 147gms |
| AU Distributor | Liberty Health | |
| Warranty | 5 Years | |



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: LSCLT_Classic-Tunable_v2-0

First Prepared: 21/01/15 Eevien Tan
 Release Approval: 21/01/15 Lee Carter
 Version Control: V2.0 Page 2 | 2

Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice.
 Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE
 Visit our website at: <http://www.lrinstruments.com.au>