

Microscope Slides & Cover Slips

German Manufactured

Microscope Slides

The German manufactured MENZEL-GLÄSER represents the best quality and experience in the world glass market today. It is associated with the most complete range of products in the field of microscope slides and cover slips.

Menzel microscope slides are available in standard sizes: 1.0 mm thick in dimensions of 26 x 76 mm (ISO Norm 8037//). Menzel microscope slides are manufactured from extra white glass. Menzel offers microscope slides with cut edges for routine microscopy applications. Menzel's single or twin frosted slides have consistent and porous surfaces making labelling effortless. The standard slide in this product line has a double frosted area of approx. 20 mm.

Microscope Cover Slips

Made of transparent hydrolytic class 1 glass, Menzel cover slips are, like other Menzel glass products, of the highest quality using glass free of blisters and flaws. With the most modern technical equipment, Menzel achieves consistency in cut, thickness and dimensional accuracy. This same technology cleans and polishes the cover slips, making them suitable for use right from the pack. Thickness 0.13mm.



CODE	DESCRIPTION
MSP	MICROSCOPE SLIDES PLAIN (PK 50)
MSF	MICROSCOPE SLIDES FROSTED (PK 50)
MSC2222	COVER SLIPS 22 X 22MM (PK 100)
MSC2240	COVER SLIPS 22 X 40MM (PK 100)

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: MS_Microscope-Slides_v2-0



First Prepared: 3/12/14 Eevien Tan
Release Approval: 3/12/14 Lee Carter
Version Control: V2.0 Page 1 | 2

Microscope Slides & Cover Slips

Chinese Manufactured

Microscope Slides

Liberty's microscope slides are pre-cleaned and suitable for any medical or research laboratory tests. The Liberty microscope slides have single frosted end on which slides have a double frosted area approximately 20mm long.



Microscope Cover Slips

The cover glass is used in conjunction with a microscope slide. One function of a cover slip is to slow evaporation from the surface of the slide being examined. The cover slip also serves to flatten the object being viewed and can serve as a barrier between this object and the microscope's objective lens to keep the lens clean.

Liberty's standard series cover slips, represent the highest quality level from China. The sterile packaging guarantees cover slips can be packed in a vacuum and dry condition. Thickness 0.13mm.



CODE	DESCRIPTION
MSPC	LIBERTY MICROSCOPE SLIDES PLAIN (PK 50)
MSFC	LIBERTY MICROSCOPE SLIDES FROSTED (PK 50)
MSC2222C	LIBERTY COVER SLIPS 22 X 22MM (PK 200)

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: MS_Microscope-Slides_v2-0



First Prepared: 3/12/14 Eevien Tan
Release Approval: 3/12/14 Lee Carter
Version Control: V2.0 Page 2 | 2