

Bestneb Nebuliser Pump

12LPM Ultra Quiet Operation

BESTNEB Nebuliser is a high quality aerosol compressor with contemporary, nonmedical styling. It is appealing to the clinician and the patient. Especially in the pediatric case, it makes children's respiratory therapy less frightening and has revolutionised the industry's approach to children's home health care products. According to the latest study, the majority of nebuliser user are from the paediatric field due to the poor resistance to air pollution.

Features

- Contemporary styling
- Ultra quiet operation
- Pleasant pastel colours
- Compartment for storing accessories
- Made from high quality Polymousse
- High Quality Aerosol Compressor
- Large Air Intake Filter

Complete with

- 1 Adult Mask
- 1 Child Mask
- 1 Nebuliser Bowl
- 1 2M supply Tubing
- 1 Mouthpiece
- 1 T Piece
- 1 Reservoir Tube 15CM

Specifications

Power:	220/240 V ~ 50 HZ
Dimensions:	29 X 20 X 12 CM
Weight:	2.5 Kg
Warranty:	2 Or 5 Years
Maximum Pressure:	32 PSI
Flow:	12 LPM
Sound Level:	50 DB
Particle Size (Mmad):	1.1 Micron



Code	Description
PNB	BESTNEB NEBULISER PUMP
PNBF	SPARE FILTER PKT 5
31310	NEBULISER TUBING 2 METRE
3569	NEBEASY NEBULISER BOWL, MASK & TUBING ADULT
3570	NEBEASY NEBULISER BOWL, MASK & TUBING CHILD

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: PNB_Nebuliser-Bestneb_v3-0



First Prepared: 11/11/14 Eevien Tan
 Release Approval: 14/11/14 Lee Carter
 Version Control: V3.0 Page 1 | 2

Bestneb Nebuliser Pump 12LPM Ultra Quiet Operation

Specific Warranty Information

- This warranty does not cover accessories, filters or damage due to misuse or negligence.
- Personal Home Use: 5 Year Warranty
- Distributors, Hospitals, Medical Centres: 2 Year Warranty
- Warranty claims or service requirement enquiries should be directed to the place of purchase.

General Product Warranty 2 Year

This warranty replaces and supersedes any other enclosed or printed warranty statement or guarantee in the instructions for use or manual.

Paragon Medical Pty Ltd trading as LR Instruments warrants that the product supplied is free from faults or defects in materials or workmanship for 2 years (24 months) from the date of original purchase.

This warranty excludes defects, faults or a failure to perform as intended caused by the product not being used in accordance with instructions, accidental or intentional damage, normal wear and tear, misuse, being altered or tampered with by unauthorised persons or any other reasons allowed by Australian Consumer Law.

The benefits provided to the consumer by this warranty are in addition to other rights and remedies available to the consumer under the law.

This warranty is given by:

Paragon Medical Pty Ltd trading as LR Instruments
Unit 7, 100 New Street, Ringwood Victoria 3134
Phone: 1800 337 048
Email: info@lrinstruments.com.au

To claim this warranty, the consumer must cease using the product immediately and, in the first instance, contact the supplier from whom the item was purchased detailing the problem and arranging the return of the faulty item for inspection (including testing if applicable) or to LR Instruments for a remedy under this warranty. Proof of purchase must accompany any claim on this warranty.

If after inspection and/or testing the product is found to be faulty or defective in materials or workmanship, the supplier or LR Instruments will at its discretion repair or replace the product. If the item is no longer available, LR Instruments reserves the right to replace it with the nearest corresponding item of no lesser value or arrange a refund through the supplier from which the item was originally purchased.

The costs of returning the faulty or defective item to the supplier or LR Instruments must be borne by the consumer.

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure.

2Year, V2.1

Warranty_General-Products_2Year_PM(LRI)_v2-1



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: PNB_Nebuliser-Bestneb_v3-0

First Prepared: 11/11/14 Eevien Tan
Release Approval: 14/11/14 Lee Carter
Version Control: V3.0 Page 2 | 2