

Liberty Nebuliser Pump

The **Liberty**[®] Nebuliser Pump is the product of 20 years of extensive research and development to provide the patient with the most cost effective, versatile and reliable pump on the market.

The **Liberty**[®] Nebuliser Pump has been designed with high performance as a priority. From the air intake filter through to the nebuliser bowl, nothing has been left to chance. The motor and precision pump assembly allows continuous operation without the annoying short-term usage times or cool down periods associated with similarly priced pumps. Use it whenever you want, how often you want and for how long you want.

Having a compact and contemporary design, the **Liberty**[®] Nebuliser Pump is a quality product suitable for use with a wide variety of medications recommended by your doctor for the treatment and/or relief of respiratory complaints.

Aerosol therapy is one of the methods of drug delivery for the treatment of chronic and acute pulmonary disease. This form of treatment is increasingly favoured by both family doctors and specialists since it achieves maximum drug efficiency directly in the respiratory tract thus minimising possible side-effects.

The **Liberty**[®] aerosol therapy system incorporates an oil free compressor. Combined with the performance of the **RAPIDFLAEM2** jet nebuliser (2.7-3µm) and a full set of accessories including mask, mouth piece and nasal prongs, the **Liberty**[®] Nebuliser Pump is a high quality system, making it suitable for basic care and seasonal allergic pathologies.

The **Liberty**[®] Nebuliser Pump has a patented anti-vibration system offering a silent and vibration-free operation. The case is moulded from thermoplastic resins, has robust features, a handy carry handle and a large accessory compartment which is hygienically separated from the cord storage compartment.

Features

- Continuously rate for constant use
- Contemporary compact styling
- Super quiet operation
- Compartment for storage of accessories

Complete with

- 1 Adult Mask
- 1 Child Mask
- 1 Nebuliser Bowl
- 1 Supply Tubing
- 1 Mouthpiece



Code	Description
PNL	LIBERTY NEBULISER PUMP
PNLF	SPARE FILTER (PKT 5)
31310	NEBULISER TUBING 2 METRE
PNLSA	NEBULISER BOWL, MASK & TUBING SET ADULT
PNLSC	NEBULISER BOWL, MASK & TUBING SET CHILD

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: PNL_Nebuliser-Liberty_v2-0



First Prepared: 11/11/14 Eevien Tan
Release Approval: 14/11/14 Lee Carter
Version Control: V2.0 Page 1 | 2

Liberty Nebuliser Pump



Specifications

Particle Size (Mmad):	2.7 - 3.0 Micron
Nebuliser Capacity:	7 MI
Power Supply:	230 V ~ 50 HZ / 130 VA
Maximum Pressure:	22 PSI
Flow:	9 LPM
Sound Level:	54 DB
Dimensions:	17 X 31 X 10 CM
Weight:	2.1 KG
Warranty:	2 Years

Useful Information

- Nebuliser bowls are disposable and need to be replaced regularly. Replace your jet nebuliser bowl every 6 months and always keep a spare in case of breakage or loss.
- Bowls should also be rinsed after each use, washed in warm water daily and allowed to dry naturally and stored in an airtight container.
- Once a week the mask, mouthpiece and bowl should be soaked in a disinfectant such as Nebuclean.

High Quality Maintenance Free Compressor
High Efficiency Filter
Approvals: Tga Artg No. 92850 Class L1a



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: PNL_Nebuliser-Liberty_v2-0

First Prepared: 11/11/14 Eevien Tan
Release Approval: 14/11/14 Lee Carter
Version Control: V2.0 Page 2 | 2

Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice.
Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE
Visit our website at: <http://www.lrinstruments.com.au>