

Electric Plaster Cutting Saw

125W Input Power and Alloy Blade



Features

- Swing: Frequency: 133Hz
- Range: ≥4 degrees
- Power Supply: 220V AC, 50Hz
- Input Power: 125W
- Blade Diameter: 62MM
- Blade Type: Oscillating (back and forth motion)
- Power Cord Length: 2M
- Low noise total optimized mechanical structure and parts design <85dB
- Alloy blade is suitable for constant use
- Function of parenchyma protection is safe and reliable
- Variable speed control
- Suitable for cutting macromolecule plastic material

Complete With

- 1 Liberty® Electric Plaster Cutting Saw
- 2 Alloy Cutting Blades 62MM
- 1 Allen Key
- 1 Lockable Attache Style Carry Case



Code	Description
PSEL	LIBERTY PLASTER SAW ELECTRIC
PSELB	LIBERTY PLASTER SAW SPARE BLADE

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: PSEL_Plaster-Saw-Electric_v3-1



First Prepared: 4/12/14
 Release Approval: 4/12/14
 Version Control: V3.1
 Eeven Tan
 Lee Carter
 Page 1 | 1