

Babinsky Telescopic Reflex Hammer

Liberty



The Liberty Babinsky Telescopic Reflex Hammer has been designed with stainless steel construction and a rubber impact ring head that can be mounted parallel to the shaft or at a right angle.

Features

- Stainless steel construction
- Rubber impact ring for testing reflex reaction
- Head can be mounted at right angle (illustrated) or parallel to the shaft
- Equipped with a pointed end piece for testing nervous system response
- Unique telescopic design makes a compact unit



Code	Description
RHBT	LIBERTY REFLEX HAMMER BABINSKY TELESCOPIC

Prepared by:	Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details:	P (03) 9870 7066 F (03) 9870 1638 E info@lrinstruments.com.au
File Name:	RHBT_Reflex-Hammers_Babinsky_v2-0



First Prepared:	19/11/14	Eevien Tan
Release Approval:	19/11/14	Lee Carter
Version Control:	V2.0	Page 1 1

Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice.
Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE
Visit our website at: <http://www.lrinstruments.com.au>