

Reflex Hammer

Percuflex

The Percuflex Reflex Hammer is specially designed with a tunable head for reflex testing and percussion.

Features

- **Handle:** The ergonomic chrome-plated brass handle is precisely balanced for increased control of force in the percussion.
- **Head:** The balance-weighted head is outfitted with a large and small mallet constructed of soft TPR for precise percussions.
- **Retractable Brush:** The built-in retractable brush at the base of the handle can be easily accessed to elicit cutaneous reflexes.
- **Integrated Pointed Tip:** The recessed chrome-plated brass head is equipped with a screw-in pointed tip for pinprick cutaneous reflex testing.



Code	Description
RHP	REFLEX HAMMER PERCUFLEX 75G BLUE

Prepared by:	Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details:	P (03) 9870 7066 F (03) 9870 1638 E info@lrinstruments.com.au
File Name:	RHP_Reflex-Hammers_Percuflex_v2-0



First Prepared:	19/11/14	Eevien Tan
Release Approval:	19/11/14	Lee Carter
Version Control:	V2.0	Page 1 1