

Electronic Blood Pressure Monitor

DS-125D With Digital Display

Features

- Measures systolic and diastolic pressure and pulse
- Easy to read digital display
- Automatic deflation
- Automatic off
- Self adjusting D-Bar cuff for one hand application

Specifications

Measurement Range	Cuff Pressure: 0 to 300mmHg Pulse: 30 to 150 pulse/minute
Accuracy	Pressure: ± 3 mm Hg Pulse: $\pm 5\%$ of reading
Temperature	+10°C to +40°C
Humidity	85% Rh or below
Storage Temperature	-30°C to +65°C
Power Supply	9 Volt Battery (006P or 6F22)
Inflation	Squeeze Bulb
Deflation	Pre-settable automatic air release system
Exhaust	Manual Valve
Cuff	Standard Adult Size
Weight	Approx. 190 g (including battery)
Size	158mm(W) x 86mm(D) x 29.6mm(H)



Product Code SDDS125
Description ELECTRONIC SPHYG WITH DIGITAL DISPLAY

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: SDDS125_Electronic-Sphyg_v3-1



First Prepared: 15/01/15 Eevien Tan
Release Approval: 15/01/15 Lee Carter
Version Control: V3.1 Page 1 | 1