



# MERCURIAL SPHYGMOMANOMETER DESK MODEL

This desk mounted mercurial sphygmomanometer is provided with a mercury locking device. It is also provided with adult velcro cuff and blood pressure bulb.

Some key features of this desk mounted mercurial sphygmomanometer are as follows:

### Features

- Rugged die cast housing
- Graduated glass tube, features 4 mm inside bore for unrestricted mercury flow
- High visible red graduations
- Locking mercury reservoir for maintenance and transportation
- Nylon, velcro cuff, bladder, bulb and valves included
- Nickel plated air release valve



CODE	DESCRIPTION
SDM	DESK MODEL MERCURIAL SPHYG
SDMLF	DESK MODEL MERCURIAL SPHYG LATEX FREE
SDMAS	DESK MODEL MERCURIAL SPHYG WITH STETHOSCOPE

