

# Surgical Skin Marker

## Sterile



### High Quality! Low Cost

The Surgical Skin Marker features a versatile tip for precise line marking, or when held at an angle a broad line marking, to allow demarcation of selected sites, prior to, or during a surgical procedure. The marker's gentian violet ink is non toxic, non tattooing and surgical prep resistant. Markings can be removed from the skin with isopropyl alcohol. Included is a 15cm, flexible ruler calibrated in both centimetres and inches.


Widely used in general surgery marking, radiotherapy and B-type ultra-sonography positioning identification and also used in medical teaching.

### Features


- The marker's gentian violet ink is recognised by CE and has been tested by an independent laboratory
- The ink is non toxic, non tattooing and surgical prep resistant
- Medium (regular) tip
- Markings can be removed from the skin with isopropyl alcohol
- Packaged with a 15cm flexible, high quality, laminated ruler both in centimetres and inches
- Meets the single use requirement to avoid infection



Packaging Sample



SMPARMT Regular is a Medium Tip (with a purple lid)



Medium (Regular) Tip

Code	Description
SMPARMT	LIBERTY SKIN MARKER & RULER, STERILE MEDIUM TIP (BX 25)

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134  
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au  
 File Name: SMPAR\_Surgical-Skin-Marker\_v3-0



First Prepared: 3/12/14 Eevien Tan  
 Release Approval: 3/12/14 Lee Carter  
 Version Control: V3.0 Page 1 | 1