

Suction Pump Accessories

Liberty



Features

- Multi-Voltage
- Fully Portable
- Oil-free lubrication pump
- Square negative pressure dial
- Built-in bottle cavity
- Low noise
- Small size (H22CM x L29CM x D20CM)
- Easy to carry
- Fully portable unit incorporating three modes of electricity: 240V AC, 12V DC and battery
- Suitable for vehicles, clinic, home and outdoor places

Technical Information

| | |
|------------------------------|-------------------------|
| Voltage: | AC 220V/240V 50HZ 12VDC |
| Power: | 90 VA |
| Max Negative Pressure Value: | ≥0.075 MPa |
| Pumping Rate: | ≥15 Litres/Minute |
| Reservoir Volume: | 1,000 mL |
| Noise: | ≤65 Db (A) |

NOTE: The canister and lid supplied with the unit are not autoclavable. They must be washed in sterilisation solution.



The 7E-D is a MULTI-VOLTAGE unit and FULLY PORTABLE

| Code | Description |
|------|--|
| SPBB | LIBERTY 7E-D SUCTION PUMP 240V/12V/BATTERY 15L |

Spares/Accessories

| | |
|---------|----------------------|
| SPBF | FILTER |
| SPBC12 | SUCTION CANNULA 12FG |
| SPBSB1L | SPARE BOTTLE 1L |

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: SPBB_Suction-Pump_v3-1



First Prepared: 19/11/14 Eevien Tan
 Release Approval: 19/11/14 Lee Carter
 Version Control: V3.1 Page 1 | 1