

7E-A Portable Suction Pump

18 litre



The Liberty 7E-A Portable Suction Unit, is a new design, oil-free lubrication pump and based on extensive experience and development trends from domestic and overseas products.

The Liberty 7E-A Suction Pump is specially designed for extracting thick liquid such as blood, phlegm, etc. Not only can this unit be used as an aspirator in the operating room or emergency room, but also as a health care product for the home.

Features

- Portable, compact small size
- Light-weight, easy to carry
- High negative pressure
- High flow
- Oil-free lubrication pump
- Square negative pressure dial
- Built-in bottle cavity
- Low noise



Technical Information

Voltage:	AC 220/240V 50HZ
Power:	≤90 VA
Max Negative Pressure Value:	≥0.075 MPa
Pumping Rate:	≥18 Litres/Minute
Reservoir Volume:	1000 mL, 1 Piece
Noise:	≤65 Db (A)
Canister:	1 Litre

NOTE: The canister and lid supplied with the unit are not autoclavable. They must be washed in sterilisation solution.



Code	Description
SPB	LIBERTY 7E-A SUCTION PUMP 240V 18L

Spares/Accessories

SPBF	FILTER
SPBC12	SUCTION CANNULA 12FG
SPBSB1L	SPARE BOTTLE 1L

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: SPB_Suction-Pump_v3-1



First Prepared:	19/11/14	Eevien Tan
Release Approval:	19/11/14	Lee Carter
Version Control:	V3.1	Page 1 1