

Portable Suction Pump

Zeiner Rescue, Multi-Voltage 25 Litre

Battery Operated Aspirator for First Aid Emergency

Features

- Ideal for ambulance, first-aid environments and minor surgery
- Internal Ni-MH battery pack, 45 minutes run time
- Car or ambulance cable and universal voltage (100 – 240VAC) adapter and internal rechargeable battery pack
- Reliable, maintenance free and silent 12VDC double piston MEDEL pump
- Shock-resistant, easy-to-carry and easy-to-clean chassis
- Suction adjustment and vacuum gauge
- 800ml disposable canister with overflow valve
- Antibacterial filter

Technical Specification

Size:	240 x 125 x 216 mm (525 x 140 x 216 mm with double canister brackets)
Weight:	2,8 Kg
Current:	3,5 A
Max vacuum:	-82 kPa; -610 mmHg
Max flow:	25 l/min
Battery	3300 mAh (Ni-MH battery cells)
Warranty:	2 years



Code	Description
SPRESCUE	ZEINER RESCUE SUCTION PUMP UNIT 240V/12V/BATTERY 25LITRE (SEE OVER FOR SPARE PARTS LISTING)

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: SPRESCUE_Rescue-Suction-Pump_v2-0



First Prepared: 19/11/14 Eevien Tan
Release Approval: 19/11/14 Lee Carter
Version Control: V2.0 Page 1 | 2

Spares Information

For SPRESCUE Suction Pump

Autoclavable Canister for SPRESCUE

- Zeiner Rescue Suction Pump spare canister with lid
- 1200ml Autoclavable



SPRESCUEB
1200ml Autoclavable Canister

Liberty Disposable Canisters for SPRESCUE

- Single use, graduated, 4 ports, over-flow valve
- 850ml and 1200ml disposable canisters



SPCDE0101
850ml Disposable Canister

SPCDE0102
1200ml Disposable Canister

Yankauer Disposable Suction Tube

- Adult 27cm
- Sterile to avoid risk of infection - individually sealed
- Available with vacuum control
- Crown tip with vent
- Made of K-resin



LTSYDC0501
Yankauer Disposable Suction Tube Adult 27cm

Suction Pump Connecting Tubing 3/16

- 3/16" tube with connector for suckers
- Single use only - disposable
- Sterile to avoid risk of infection - individually sealed
- Blue cone connectors



SPCTC0601 Connection Tubing 3/16" ID 1.8 Metre Length
SPCTC0602 Connection Tubing 3/16" ID 3.6 Metre Length
SPCTC0603 Connection Tubing 4/16" ID 1.8 Metre Length

Code	Description
SPRESCUEB	ZEINER AUTOCLAVABLE CANISTER 1200ML FOR SPRESCUE
SPCDE0101	LIBERTY DISPOSABLE SUCTION PUMP CANISTER 850ML TO FIT SPRESCUE
SPCDE0102	LIBERTY DISPOSABLE SUCTION PUMP CANISTER 1200ML TO FIT SPRESCUE
LTSYDC0501	LIBERTY YANKAUER DISPOSABLE SUCTION TUBE ADULT 27CM
SPCTC0601	LIBERTY CONNECTION TUBING 3/16" ID 1.8 METRE LENGTH (BOX 50)
SPCTC0602	LIBERTY CONNECTION TUBING 3/16" ID 3.6 METRE LENGTH (BOX 50)
SPCTC0603	LIBERTY CONNECTION TUBING 4/16" ID 1.8 METRE LENGTH (BOX 50)
SPRESCUTBAT	BATTERY FOR SPRESCUE (NOT PICTURED)
SPRESCUEF	FILTER FOR SPRESCUE (NOT PICTURED)
SPRESCUET	TRANSFORMER FOR SPRESCUE (NOT PICTURED)



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: SPRESCUE_Rescue-Suction-Pump_v2-0

First Prepared: 4/12/14 Eevien Tan
Release Approval: 4/12/14 Lee Carter
Version Control: V2.0 Page 2 | 2

Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice.
Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE
Visit our website at: <http://www.lrinstruments.com.au>