

VacMaxi 240V

Suction Pump Powerful Suction

Smooth hygienic finish, water resistant, easy to clean and maintain.

The VAC-series pump is Apex Medical's most recent generation of portable suction units and has been designed to combine function with aesthetic appeal and meet the needs of today's hygiene and infection control requirements.

The VacMaxi suction unit is a powerful, fast and effective aspirator which is ideal for use in high vacuum, high flow applications.

Specifications

Electrical Rating:	AC 220~240V, 50/60Hz
Power Consumption:	100W
Canister Size:	1800 ML Standard
Maximum Pressure:	600 mmHg (adjustable)
Maximum Flow Rate:	46 L.P.M.
Net Weight:	4.4 kg
Dimensions:	352(L) x 206(W) x 192(H) mm
Sound Level:	< 60 dB
Warranty:	2years

Conversion: 1 millilitre = 1 cubic centimetre

PEX
VacMaxi 0197

Approval:
EN 60601-1
EN 60601-1-2
EN/ISO 10079-1

Integrated Canister Holder.

Portable design with ergonomic easy-to-carry handle.

Pressure regulator and easy-to-read analogue pressure gauge.

Detachable power cord allows you to quickly change to another cord to match local electrical outlet requirement.

Code	Description
SPVACMAXI	VACMAXI 240V SUCTION PUMP 46 LPM C/W 1800ML CANISTER
SPVACSB1800	SUCTION PUMP CANISTER 1800 ML
SPCDE0101	SUCTION PUMP CANISTER DISPOSABLE 850 ML
SPCDE0302	SUCTION PUMP CANISTER DISPOSABLE 1200 ML
SPCTC0601	SUCTION PUMP TUBING 3/16 WITH CONNECTOR 1.8M STERILE
SPCTC0602	SUCTION PUMP TUBING 3/16 WITH CONNECTOR 3.6M STERILE
SPCTC0603	SUCTION PUMP TUBING 4/16 WITH CONNECTOR 1.8M STERILE
SPBF	SUCTION PUMP FILTER

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134

Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au

File Name: SPVACMAXI_Suction-Pump_v3-0



First Prepared: 19/11/14 Eevien Tan
Release Approval: 19/11/14 Lee Carter
Version Control: V3.0 Page 1 | 1