

# VacPro 240V

## Suction Pump Lightweight

Smooth hygienic finish, water resistant, easy to clean and maintain.

The VAC-series pump is Apex Medical's most recent generation of portable suction units and has been designed to combine function with aesthetic appeal and meet the needs of today's hygiene and infection control requirements.

The VacPro is durable, lightweight and versatile in the provision of general suction usage.

### Specifications

Electrical Rating:	AC 220~240V, 50/60Hz
Power Consumption:	90W
Canister Size:	800 ML Standard or 800 ML plus 1800 ML option
Maximum Pressure:	600 mmHg (adjustable)
Maximum Flow Rate:	24 L.P.M.
Net Weight:	2.8 kg
Dimensions:	352(L) x 206(W) x 192(H) mm
Sound Level:	< 55 dB
Warranty:	2 years



Code	Description
SPVACPRO	VACPRO 240V SUCTION PUMP LIGHTWEIGHT C/W 800ML CANISTER
SPVACPRO1800	VACPRO 240V SUCTION PUMP LIGHTWEIGHT C/W 800ML & 1800ML CANISTERS
SPVACSB1800	SUCTION PUMP CANISTER 1800 ML
SPCDE0101	SUCTION PUMP CANISTER DISPOSABLE 850 ML
SPCDE0302	SUCTION PUMP CANISTER DISPOSABLE 1200 ML
SPCTC0601	SUCTION PUMP TUBING 3/16 WITH CONNECTOR 1.8M STERILE
SPCTC0602	SUCTION PUMP TUBING 3/16 WITH CONNECTOR 3.6M STERILE
SPCTC0603	SUCTION PUMP TUBING 4/16 WITH CONNECTOR 1.8M STERILE
SPBF	SUCTION PUMP FILTER

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134

Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au

File Name: SPVACPRO\_Suction-Pump-Lightweight\_v3-0



First Prepared: 19/11/14 Eevien Tan  
Release Approval: 19/11/14 Lee Carter  
Version Control: V3.0 Page 1 | 1