

# Aneroid Sphygmomanometer

## Two-Hand

### Features

- Black enamel 300 mmHg no stop pin manometer
- Velcro cuff, complete with bladder, bulb and valve
- Chrome plated deflation valve
- Comes as standard with adult cuff and bladder 2-tubes
- Made in Japan
- Complete with leatherette zipped case
- Cuffs also available in infant, child, large adult and obese as separate purchases
- Two-handed use, suitable for both right and left hands
- Range 0-300 mmHg



Comes with Adult Cuff as Standard  
Available in Hang Sell Packaging.  
To order, please add a CS to the end of the Product Code below.

Product Code	Description
SSAVC	HAND HELD ANEROID SPHYG (NAVY BLUE)
SSAVCLF	HAND HELD ANEROID SPHYG LATEX FREE (NAVY BLUE)
SSAVCBK	HAND HELD ANEROID SPHYG BLACK
SSAVCBKLF	HAND HELD ANEROID SPHYG BLACK LATEX FREE
SSAVCM	HAND HELD ANEROID SPHYG MAGENTA
SSAVCMLF	HAND HELD ANEROID SPHYG MAGENTA LATEX FREE
SSAVCNB	HAND HELD ANEROID SPHYG NAVY BLUE
SSAVCNBLF	HAND HELD ANEROID SPHYG NAVY BLUE LATEX FREE
SSAVCPU	HAND HELD ANEROID SPHYG PURPLE
SSAVCPULF	HAND HELD ANEROID SPHYG PURPLE LATEX FREE
SSAVCRB	HAND HELD ANEROID SPHYG ROYAL BLUE
SSAVCRBLF	HAND HELD ANEROID SPHYG ROYAL BLUE LATEX FREE
SSAVCT	HAND HELD ANEROID SPHYG TEAL
SSAVCTLF	HAND HELD ANEROID SPHYG TEAL LATEX FREE



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107  
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134  
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au  
File Name: SSAVC\_Sphyg-Aneroid-Two-Hand\_v3-1

First Prepared: 15/01/15 Eevien Tan  
Release Approval: 15/01/15 Lee Carter  
Version Control: V3.1 Page 1 | 2

Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice.  
Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE  
Visit our website at: <http://www.lrinstruments.com.au>