

Scalpel Blades / Handles


Sterile



These Paramount Surgical blades are available in Stainless Steel.

Each blade is protected by VCI paper in a peel-open foil sachet. It is sterilised by Gamma Radiation of 2.5 M Rads.

Scalpel Blades



Product Code	Description	Size
SSB10	NO. 10 SCALPEL BLADE	10
SSB10A	NO. 10A SCALPEL BLADE	10A
SSB11	NO. 11 SCALPEL BLADE	11
SSB12	NO. 12 SCALPEL BLADE	12
SSB15	NO. 15 SCALPEL BLADE	15
SSB20	NO. 20 SCALPEL BLADE	20
SSB21	NO. 21 SCALPEL BLADE	21
SSB22	NO. 22 SCALPEL BLADE	22
SSB23	NO. 23 SCALPEL BLADE	23
SSB24	NO. 24 SCALPEL BLADE	24
SSB25	NO. 25 SCALPEL BLADE	25 (NOT PICTURED)

Scalpel Handles



Product Code	Description	Fits
Superior Theatre Instruments		
HTH-095-003-PMK	NO. 3 SCALPEL HANDLE	BLADES 10-15
HTH-096-004-PMK	NO. 4 SCALPEL HANDLE	BLADES 20-25
HTH-097-003-PMK	NO. 7 SCALPEL HANDLE	BLADES 10-15
Theatre Instruments		
AD010/03	NO. 3 SCALPEL HANDLE	BLADES 10-15
AD010/04	NO. 4 SCALPEL HANDLE	BLADES 20-25
AD025/07	NO. 7 SCALPEL HANDLE	BLADES 10-15
Baron		
HSHB	STAINLESS STEEL HANDLE	BLADES 10-15
HSHBA	ALLOY HANDLE	BLADES 10-15
Liberty® Superior		
LSHN3	NO. 3 SCALPEL HANDLE	BLADES 10-15
LSHN4	NO. 4 SCALPEL HANDLE	BLADES 20-25
L.R.I. Ward		
SHN3	NO. 3 SCALPEL HANDLE	BLADES 10-15
SHN4	NO. 4 SCALPEL HANDLE	BLADES 20-25
SHN7	NO. 7 SCALPEL HANDLE	BLADES 10-15
SHN9	NO. 9 SCALPEL HANDLE	BLADES 10-15



Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: SSB_SSH_Scalpel-Blades-Handles_v3-1

First Prepared: 19/11/14 Eevien Tan
 Release Approval: 19/11/14 Lee Carter
 Version Control: V3.1 Page 1 | 1

Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice. Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE
 Visit our website at: <http://www.lrinstruments.com.au>