

Disposable Biopsy Punches

Sterile - Ready To Use



Standard disposable sterile scalpels offer convenience in a number of user environments be it on the ward, in emergency or trauma units or even out in the field with the paramedics.

These blades consist of pre-assembled stainless steel surgical blades mounted on an ABS plastic handle with fluted finger grips. The blades are safely protected by LDPE guard.

These scalpels are individually packed in a peel-open paper/poly pouch packed in a box of 10 sterile, disposable scalpels.

Package: 10 sterile scalpels per dispenser box



10 11 15 20 21 22 23



Advantages of using the surgical disposable scalpel

Using a sterile disposable scalpel can reduce or even prevent the incidences of:

- Injuries during the removal of blades from re-usable handles
- Occupational illness arising from cross infection (particularly of highly infectious and life threatening diseases such as HIV/ Aids, Hepatitis B, C) caused through injuries sustained during the removal of blades from re-usable handles
- These scalpels are meant for one time use only. The sterility of a freshly opened disposable scalpel is assured.
- After using the disposable scalpel(s), the contaminated blades can be easily covered with the protective caps or guards, then safely disposed of in the Operation Theater/ Room Disposal Container. There is no risk of puncture, due to this convenient and easy way of disposing the sharp blades.
- Being a ready-to-use "sterile" product is useful for immediate surgery.

| Product Code | Description |
|--------------|---|
| SSDC10 | DISPOSABLE SCALPEL COMPLETE, NO. 10 (BOX 10) |
| SSDC11 | DISPOSABLE SCALPEL COMPLETE, NO. 11 (BOX 10) |
| SSDC12 | DISPOSABLE SCALPEL COMPLETE, NO. 12 (BOX 10) [NOT PICTURED] |
| SSDC15 | DISPOSABLE SCALPEL COMPLETE, NO. 15 (BOX 10) |
| SSDC20 | DISPOSABLE SCALPEL COMPLETE, NO. 20 (BOX 10) |
| SSDC21 | DISPOSABLE SCALPEL COMPLETE, NO. 21 (BOX 10) |
| SSDC22 | DISPOSABLE SCALPEL COMPLETE, NO. 22 (BOX 10) |
| SSDC23 | DISPOSABLE SCALPEL COMPLETE, NO. 23 (BOX 10) |
| SSDC24 | DISPOSABLE SCALPEL COMPLETE, NO. 24 (BOX 10) [NOT PICTURED] |

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
 Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
 File Name: SSDC_Scalpels-Disposable_v3-1



First Prepared: 19/11/14 Eevien Tan
 Release Approval: 19/11/14 Lee Carter
 Version Control: V3.1 Page 1 | 1