

Aneroid Sphygmomanometer With Trolley & Basket



Sphygmomanometer Features

- Large, easy to read 150mm dial
- 300 mm Hg no pin stop mechanism
- Rear bracket storage compartment area
- Full 110° rotation for easy viewing
- Nylon velcro cuff, bladder, bulb and valve included
- Nickel plated air release valve
- 8 foot (when extended) of coil tubing

Trolley, Basket Features

- Basket is a powder coated hard wire mesh
- Stainless steel pole
- Easy flow, high quality castors
- Broad base for extra stability
- Base options in light weight, medium weight and heavy weight



Powder coated hard wire mesh basket



Large, easy to read 150 mm dial



Code	Description
STMALW	ANEROID SPHYG WITH TROLLEY & BASKET LIGHT WEIGHT BASE
STMALWLF	ANEROID SPHYG WITH TROLLEY & BASKET LIGHT WEIGHT BASE LATEX FREE
STMAMW	ANEROID SPHYG WITH TROLLEY & BASKET MEDIUM WEIGHT BASE
STMAMWLF	ANEROID SPHYG WITH TROLLEY & BASKET MEDIUM WEIGHT BASE LATEX FREE
STMAHW	ANEROID SPHYG WITH TROLLEY & BASKET HEAVY WEIGHT BASE
STMAHWLF	ANEROID SPHYG WITH TROLLEY & BASKET HEAVY WEIGHT BASE LATEX FREE
SCT	COILED TUBING FOR SPHYG (2 METRES) (NOT PICTURED)

Prepared by: Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107
Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details: P (03) 9870 7066 | F (03) 9870 1638 | E info@lrinstruments.com.au
File Name: STMA_Aneroid-Sphyg_v3-2



First Prepared: 15/01/15 Eevien Tan
Release Approval: 15/01/15 Lee Carter
Version Control: V3.2 Page 1 | 1