



TROLLEY & WALL MODEL MERCURIAL SPHYGMOMANOMETER

Features

- Rugged die cast housing
- Graduated glass tube, features 4 mm inside bore for unrestricted mercury flow
- High visible red graduations
- Locking mercury reservoir for maintenance and transportation
- Nylon, velcro cuff, bladder, bulb and valves included
- 8 foot (when extended) of coiled tubing
- Nickel plated air release valve
- 4 or 5 castor stand (on trolley model)
- Telescopic pole (on trolley model)
- Basket for cuff storage (on trolley model)



The Wall Model also features a full 110° swivel for enhanced visibility and view comfort



CODE	DESCRIPTION
STMM4	TROLLEY MODEL MERCURIAL 4 LEG
STMM4LF	TROLLEY MODEL MERCURIAL 4 LEG - LATEX FREE
STMM5	TROLLEY MODEL MERCURIAL 5 LEG
STMM5LF	TROLLEY MODEL MERCURIAL 5 LEG - LATEX FREE
SWMM	WALL MOUNTABLE MERCURIAL SPHYG
SWMLF	WALL MOUNTABLE MERCURIAL SPHYG - LATEX FREE

