

# Pediatric Urine Collectors



In order to obtain suitable specimens from infants and small children, paediatric urine collectors which are attached to the genitalia, are available from Liberty.

## Features

- Foam ring at opening for patient comfort
- 100cc graduated
- Adhesive surface for security and and sealing container
- Suitable for male and female patients
- Available in box of 100 collectors.



Code	Description
UCP100	URINE COLLECTORS PAEDIATRIC (BOX 100)

Prepared by:	Paragon Medical Pty Ltd trading as LR Instruments ABN 76 133 072 107 Unit 7/100 New Street (PO Box 458) Ringwood VIC 3134
Contact Details:	P (03) 9870 7066   F (03) 9870 1638   E info@lrinstruments.com.au
File Name:	UCP_Pediatric-Urine-Collector_v2-0



First Prepared:	13/01/15	Eevien Tan
Release Approval:	13/01/15	Lee Carter
Version Control:	V2.0	Page 1   1

Disclaimer: All images are included for reference only. Features and specifications are subject to change without prior notice.  
Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes. E & OE  
Visit our website at: <http://www.lrinstruments.com.au>