



LIGHTWEIGHT MANUAL WHEELCHAIR

ALUMINIUM FRAME - ONLY WEIGHS 14KG

- Lightweight aluminium frame
- Only weighs 14kg when ready for transport
- Tipping levers for attendant
- Easy to fold for convenient transport or storage
- Rear wheels 610mm with strong plastic handrims
- Fixed arm rests and foot rests
- Dual side brakes

Technical Information

Maximum Load	100KG
Seat Width	46CM
Seat Depth	40CM
Seat Height	40CM
Attendant Handle Height	87CM
Wheelchair Length	100CM
Wheelchair Width	65CM
Wheelchair Width Folded	24CM
Front Wheel	20CM Castor
Rear Wheel	61CM PU Wheel

**new
product**

**wheelchair
only weighs
14kg**



CODE
WCAKY868L

DESCRIPTION
WHEELCHAIR LIGHTWEIGHT ALUMINIUM

E&OE First Published: 17 March 2011
File Name: info_WCAKY868L_WHEELCHAIR-LIGHTWEIGHT

