

TRANSIT WHEELCHAIR

Simple, robust and very low cost transit chair. Effortless to push and comfortable to ride. A lightweight push/transit wheelchair featuring a folding, steel frame designed for those users unable to propel themselves.

This attendant propelled Transit wheelchair folds making it easily transported in most cars and very convenient for home care personnel.

Features:

- Lightweight design
- Removeable foot plates
- Removeable armrests
- Dual side brakes
- Standard mag rear wheel with pneumatic tyre
- Folding frame for convenient transport
- Synthetic, washable padded upholstery
- Tipping levers for the attendant
- Burgundy colour

Technical Data:

- Weight: 16kg
- Maximum Load: 100kg
- Seat Width: 46cm
- Seat Depth: 40cm



CODE	DESCRIPTION
WCT100B	TRANSIT WHEELCHAIR
WCST	SPARE TUBE FOR TRANSIT REAR WHEEL

E&OE First Published: 16 July 2009
Last Updated: 14 February 2011
File Name: info_WCT100B_TRANSIT-WHEELCHAIR