



Liberty®

LIBERTY NEBULISER PUMP

Continuous Use, No Resting

The Liberty® Nebuliser Pump is the product of 20 years of extensive research and development to provide the patient with the most cost effective, versatile and reliable pump on the market.

The Liberty® Nebuliser Pump has been designed with high performance as a priority. From the air intake filter through to the nebuliser bowl, nothing has been left to chance. The motor and precision pump assembly allows continuous operation without the annoying short-term usage times or cool down periods associated with similarly priced pumps. Use it whenever you want, how often you want and for how long you want.

Complete with

- 1 Adult Mask
- 1 Child Mask
- 1 Nebuliser Bowl
- 1 Supply Tubing
- 1 Mouthpiece

Features

- Continuously rate for constant use
- Contemporary compact styling
- Super quiet operation
- Compartment for storage of accessories

Having a compact and contemporary design, the Liberty® Nebuliser Pump is a quality product suitable for use with a wide variety of

medications recommended by your doctor for the treatment and/or relief of respiratory complaints.

Aerosol therapy is one of the methods of drug delivery for the treatment of chronic and acute pulmonary disease. This form of treatment is increasingly favoured by both family doctors and specialists since it achieves maximum drug efficiency directly in the respiratory tract thus minimising possible side-effects.

The Liberty® aerosol therapy system incorporates an oil free compressor. Combined with the performance of the RAPIDFLAEM2 jet nebuliser (2.7-3µm) and a full set of accessories including mask, mouth piece and nasal prongs, the Liberty® Nebuliser Pump is a high quality system, making it suitable for basic care and seasonal allergic pathologies.

The Liberty® Nebuliser Pump has a patented anti-vibration system offering a silent and vibration-free operation. The case is moulded from thermoplastic resins, has robust features, a handy carry handle and a large accessory compartment which is hygienically separated from the cord storage compartment.

CODE	DESCRIPTION
PNL	LIBERTY NEBULISER PUMP
31310	NEBULISER TUBING 2 METRE
PNLF	SPARE FILTER (PKT 5)
PNLMP	SPARE MOUTH PIECE (#10 IN INSTRUCTION MANUAL)
PNLMS	SPARE MEDICINE SAVER (#7 IN INSTRUCTION MANUAL)
PNLSA	NEBULISER BOWL, MASK AND TUBING ADULT
PNLSC	NEBULISER BOWL, MASK AND TUBING CHILD

All images are included for reference only. Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes.



BESTNEB NEBULISER PUMP

Use for 25 minutes, then rest!

BESTNEB nebuliser is a high quality aerosol compressor with contemporary, non medical styling. It is appealing to the clinician and the patient. Especially in the paediatric case, it makes children's respiratory therapy less frightening and has revolutionised the industry's approach to children's home health care products. According to the latest study, the majority of nebuliser user are from the paediatric field due to the poor resistance to air pollution.

Complete with

- 1 Adult Mask
- 1 Child Mask
- 1 Nebuliser Bowl
- 1 2M supply Tubing
- 1 Mouthpiece
- 1 T Piece
- 1 Reservoir Tube 15CM

Features

- Contemporary styling
- Ultra quiet operation
- Pleasant pastel colours
- Compartment for storing accessories

Specifications

POWER:	220/240 V ~ 50 HZ
DIMENSIONS:	29 X 20 X 12 CM
WEIGHT:	2.5 KG
WARRANTY:	2 OR 5 YEARS
MAXIMUM PRESSURE:	32 PSI
FLOW:	12 LPM
SOUND LEVEL:	50 DB
PARTICLE SIZE:	.5 - 5 MICRON
HIGH QUALITY AEROSOL COMPRESSOR	
LARGE AIR INTAKE FILTER	

CODE	DESCRIPTION
PNB	BESTNEB NEBULISER PUMP
31310	NEBULISER TUBING 2 METRE
3569	NEBEASY NEBULISER BOWL, MASK & TUBING ADULT
3570	NEBEASY NEBULISER BOWL, MASK & TUBING CHILD
PNBF	SPARE FILTER PKT 5

NEBULISER PUMP FOOT

EASY AIR Clinical Nebuliser Pump

- 100% portable and independent
- No batteries or power packs required
- Lightweight (700 gm)
- Small (111 x 190 x 306 mm)
- High capacity (9.4 Litre/min @ 50 strokes/min)
- Wide range of operating temperatures (5 - 35°C)
- Robust



CODE	DESCRIPTION
PNF	FOOT/HAND NEBULISER PUMP

All images are included for reference only. Goods supplied may vary from those illustrated due to continual improvements and/or regulation changes.