

Nebuliser Pump | Liberty



Continuous Use, No Resting

The Liberty® Nebuliser Pump is the product of 20 years of extensive research and development to provide the patient with the most cost effective, versatile and reliable pump on the market.

The Liberty® Nebuliser Pump has been designed with high performance as a priority. From the air intake filter through to the nebuliser bowl, nothing has been left to chance. The motor and precision pump assembly allows continuous operation without the annoying short-term usage times or cool down periods associated with similarly priced pumps. Use it whenever you want, as often as you want and for as long as you want.

Having a compact and contemporary design, the Liberty® Nebuliser Pump is a quality product suitable for use with a wide variety of medications recommended by your doctor for the treatment and/or relief of respiratory complaints.

Aerosol therapy is one of the methods of drug delivery for the treatment of chronic and acute pulmonary disease. This form of treatment is increasingly favoured by both family doctors and specialists since it achieves maximum drug efficiency directly in the respiratory tract thus minimising possible side-effects.

The Liberty® aerosol therapy system incorporates an oil free compressor. Combined with the performance of the RAPIDFLAEM2 jet nebuliser (2.7-3µm) and a full set of accessories including mask, mouth piece and nasal prongs, the Liberty® Nebuliser Pump is a high quality system, making it suitable for basic care and seasonal allergic pathologies.

Product Code	Description
PNL	LIBERTY NEBULISER PUMP
31310	NEBULISER TUBING 2 METRE
PNLF	SPARE FILTER (PKT 5)
PNLMS	SPARE MEDICINE SAVER (#7 IN INSTRUCTION MANUAL)
PNLSA	NEBULISER BOWL, MASK AND TUBING ADULT
PNLSC	NEBULISER BOWL, MASK AND TUBING CHILD

The Liberty® Nebuliser Pump has a patented anti-vibration system offering a silent and vibration-free operation. The case is moulded from thermoplastic resins, has robust features, a handy carry handle and a large accessory compartment which is hygienically separated from the cord storage compartment.

Features

- "Continuously" rated for constant use
- Contemporary compact styling
- Super quiet operation
- High quality maintenance free compressor
- High efficiency filter
- Made in Italy
- Compartment for storage of accessories

Complete with

- 1 Adult Mask
- 1 Child Mask
- 1 Nebuliser Bowl
- 1 Supply Tubing
- 1 Mouthpiece
- 1 Nasal Prongs

Specifications

- Particle Size (mmAD): 2.7 - 3.0 micron
- Capacity: 7 ml
- Power Supply: 230 v ~ 50 Hz / 130 vA
- Maximum Pressure: 22 psi
- Flow: 9 lpm
- Sound Level: 54 dB
- Dimensions: 17 x 31 x 10 cm
- Weight: 2.1 kg
- Warranty: 2 year manufacturers warranty

Nebuliser Pump | Bestneb

**Use for 25 minutes, then rest!**

Bestneb nebuliser is a high quality aerosol compressor with contemporary, non medical styling. It is appealing to the clinician and the patient. Especially in pediatric cases, it makes children's respiratory therapy less frightening and has revolutionised the industry's approach to children's home health care products. According to the latest study, the majority of nebuliser users are from the pediatric field due to the poor resistance to air pollution.

This unit is intended for intermittent use only (max 25 minutes at a time). If you require a "continuously" rated unit, see the Liberty Nebuliser Pump.

Features

- Contemporary styling in pleasant pastel colours
- Ultra quiet operation
- High quality aerosol compressor
- Large air intake filter
- Made in Taiwan
- Compartment for storing accessories

Comes complete with:

1 Adult Mask, 1 Child Mask, 1 Nebuliser Bowl, 1 x 2m Supply of Tubing, 1 Mouthpiece, 1 T-piece, 1 Reservoir Tube 15cm

Warranty

- 2 year warranty for professional use (hospitals, clinics, hiring etc)
- Extended 5 year warranty for home use only

Specifications

- Particle Size (mmAD): 1.1 micron
- Power Supply: 220/240 v ~ 50 Hz
- Maximum Pressure: 32 psi
- Flow: 12 lpm
- Sound Level: 55 db
- Dimensions: 29 x 20 x 12 cm
- Weight: 2.5 kg

Product Code

PNB	DESCRIPTION
31310	BESTNEB NEBULISER PUMP
3551	NEBULISER TUBING 2 METRE
3569	NEBEASY JET NEBULISER BOWL ONLY
3570	NEBEASY NEBULISER BOWL, MASK & TUBING ADULT
PNBF	NEBEASY NEBULISER BOWL, MASK & TUBING CHILD
75047	SPARE FILTER PKT 5
	SPARE T-PIECE

Aerosol Masks



Liberty[®]

- Soft, transparent, odour free vinyl mask
- Swivel connection
- 2 metre tubing
- Jet nebuliser bowl
- Graduated bowl
- Dead Space: 0.5ml
- 0.5 - 5 micron
- Flow Range: 4 - 12 lpm



Product Code
LRESA0205
LRESA0206


Description
LIBERTY JET NEBULISER BOWL/MASK/TUBING ADULT
LIBERTY JET NEBULISER BOWL/MASK/TUBING CHILD

Oxygen Masks



Liberty[®]

- Soft, transparent, odour free vinyl
- Reinforced strip to ensure the best fit and user comfort
- Swivel connection



Product Code
LRESA0105
LRESA0106

Description
LIBERTY OXYGEN MASK WITH TUBING ADULT
LIBERTY OXYGEN MASK WITH TUBING CHILD

Hi-Oxygen Concentration Masks



Liberty[®]

- Soft, transparent, odour free vinyl
- Aluminium fixer and elastic strip to ensure the best fit and user comfort
- One way valve and reservoir to increase the efficiency of oxygen concentration
- Swivel connection



Product Code
LRESA0301
LRESA0302

Description
LIBERTY HIGH OXYGEN MASK, TUBING/RESERVOIR ADULT
LIBERTY HIGH OXYGEN MASK, TUBING/RESERVOIR LARGE CHILD/STANDARD ADULT

Oxygen Bottle & Gas Cylinder Trolleys



Oxygen Bottle Trolleys

- Epoxy coated or stainless steel
- 150mm Castors



Multi Stack Gas Cylinder Trolleys

- Epoxy coated
- C Size cylinders
- 2 different sized trolleys with carrying capacity of 9 or 12 cylinders

ParagonCare

paragoncare.com.au

Oxygen Cannula

- Soft, transparent, odour free PVC to ensure the best fit and user comfort
- Soft tip, over the ear style to provide low flow oxygen
- Disposable, anti crush 2 metre tubing
- Latex Free




Product Code
LRESC0205
LRESC0206
LRESC0207
LRESC0208

Description
LIBERTY SOFT NASAL CANNULA 2M TUBING ADULT
LIBERTY SOFT NASAL CANNULA 2M TUBING CHILD
LIBERTY SOFT NASAL CANNULA 2M TUBING INFANT
LIBERTY SOFT NASAL CANNULA 2M TUBING NEONATE

Oxygen Tubing

- Disposable oxygen supply tubing, crush-resistant
- Length = 2 metres




Product Code
LRESC0401

Description
LIBERTY OXYGEN TUBING 2 METRES

Oxygen Tubing Connector

- Oxygen supply tubing connector



Product Code
LRESPN1139

Description
LIBERTY OXYGEN TUBING CONNECTOR (PACK 20)

