



BACTERIAL FILTER WITH O2 PORT

Features

- Designed for a better inlet air quality
- Designed to prevent cross contamination and protection of the operator
- Efficiency: 99% for 0.3 micron (3M's material)
- Connector 22 mm ID and 22/15 mm



CODE	DESCRIPTION
39016	DISPOSABLE BACTERIAL FILTER WITH LEUR PORT
39014	HEPA RADIAL PLEATED FILTER
70015	BACTERIAL FILTER WITH CLEAR HOUSING (NOT PICTURED)

E&OE First Published: 7 April 2009
Last Modified: 25 November 2009
File Name: info_3901X_BACTERIAL-FILTERS

