



OXYGEN MASKS

MASKS - ECONOMY SILICONE

- Material: Polycarbonate and/or silicone rubber
- Autoclavable up to 134°C
- Connections: 15mm OD for 5121, #5120
22mm ID for 5122, 5133, 5134 and 5135

CODE	DESCRIPTION
5120	ECONOMY SILICONE MASK NO. 0 INFANT
5121	ECONOMY SILICONE MASK NO. 1 LARGE INFANT
5122	ECONOMY SILICONE MASK NO. 2 CHILD
5133	ECONOMY SILICONE MASK NO. 3 LARGE CHILD
5134	ECONOMY SILICONE MASK NO. 4 ADULT
5135	ECONOMY SILICONE MASK NO. 5 LARGE ADULT



MASKS - DURABLE SILICONE

- Material: Polysulphone and silicone rubber
- Fully autoclavable up to 170°C
- Connections - 22 mm ID
- Dead Space: Adult mask less than 150 ml;
Child mask less than 95 ml

CODE	DESCRIPTION
5143	DURABLE SILICONE MASK NO. 3 LARGE CHILD
5144	DURABLE SILICONE MASK NO. 4 ADULT
5145	DURABLE SILICONE MASK NO. 5 LARGE ADULT
15-001	SILICONE SIDE RETAINER FOR MASKS (PACK 10)
15-002	SILICONE FRONT RETAINER FOR MASKS (PACK 10)
15-005	SILICONE TOP GROMMET FOR MASK (PACK 5)
15-003	SILICONE CUSHION FOR 5145 MASK
15-007	SILICONE CUSHION FOR 5144 MASK
15-009	SILICONE CUSHION FOR 5143 MASK



MASKS - RENDELL BAKER

- Transparent silicone rubber
- 22mm Connection
- Minimises dead space
- Autoclavable up to 170°C

CODE	DESCRIPTION
5150	SILICONE RENDELL BAKER MASK NO. 0 INFANT
5151	SILICONE RENDELL BAKER MASK NO. 1 CHILD SMALL
5152	SILICONE RENDELL BAKER MASK NO. 2 CHILD MEDIUM
5153	SILICONE RENDELL BAKER MASK NO. 3 CHILD LARGE



E&OE First Published: 20 March 2009
Last Modified: 24 November 2009
File Name: info_51XX_OXYGEN-MASKS

