



CPAP MASK WITH HEAD HARNESS

Features

- CPAP Full Face Mask offers continue positive airway pressure for treatment of obstructive sleep apnea, a disorder which disrupts a person's breathing during sleep
- Silicone cushion are made of soft medical grade silicone material for better seal and for patient's comfort
- Double swivel elbow increases the flexibility and comfort for the patient during sleep
- Build in 4.0M pressure port features
 - a) pressure sensing,
 - b) CO2 wash out,
 - c) oxygen supply
- Crown and Easy release strap provides extra stability, eliminated the need for a chin strap, and the easy-release code for one-step removal of headgear



Specifications

- Material conform to ISO-10993 for biocompatibility requirement
- 22M connector for tubing, and 4.0M for pressure sensing, CO2 sensing, and oxygen supply

CODE	DESCRIPTION
5427	CPAP FULL FACE MASK LARGE WITH HARNESS
5428	CPAP FULL FACE MASK MEDIUM WITH HARNESS
5429	CPAP FULL FACE MASK SMALL WITH HARNESS

E&OE First Published: 20 March 2009
Last Modified: 24 November 2009
File Name: info_542X_CPAP-WITH-HEAD-HARNESS

