



VM-2 ANIMAL MASKS

GaleMed VM-2 Animal Mask is designed for small animals for anesthesia, resuscitation, oxygen therapy and aerosol therapy. Made of clear high impact polycarbonate for shell, stainless steel hook ring, and silicone rubber for cushion, harness and 22F universal connection. 6 different sizes available for different small animals and a standard 15M or 22F connection for ventilation equipment.



Features

- Dimensions: 15mm OD for circuit connection
- Cushion ID: 33mm
- Flexibility of the diaphragm ensures air tight seal
- Shape of cone ensures low air resistance
- Transparent cone permits clear visibility of patients muzzle
- Material: polycarbonate and silicone rubber
- Very low dead space, important when masking down using circle systems
- Removeable rubber diaphragm which minimises gas leaks and permits cleaning
- Made from virtually indestructible PVC which will not degrade

CODE

0011	SILICONE HEAD HARNESS - LARGE
0013	SILICONE HEAD HARNESS - SMALL
1580-02	CUSHION FOR ANIMAL MASK - SS
1581-02	CUSHION FOR ANIMAL MASK - SM
1582-02	CUSHION FOR ANIMAL MASK - S
1583-02	CUSHION FOR ANIMAL MASK - M
1584-02	CUSHION FOR ANIMAL MASK - L
1585-02	CUSHION FOR ANIMAL MASK - LL
(CUSHIONS NOT PICTURED)	
5630	VM-2 ANIMAL MASK #0, SS 50MM/26MM ID
5631	VM-2 ANIMAL MASK #1, SM 60MM/31MM ID
5632	VM-2 ANIMAL MASK #2, S 65MM/39MM ID
5633	VM-2 ANIMAL MASK #3, M 95MM/36MM ID
5634	VM-2 ANIMAL MASK #4, L 115MM/43MM ID
5635	VM-2 ANIMAL MASK #5, LL 135MM/74MM ID

E&OE First Published: 20 March 2009
 Last Modified: 24 November 2009
 File Name: info_563X_VM2-ANIMAL-MASKS

