



CLINIPAK CORNERS


Perfect partner for Loaset Management!

A simple and effective solution to reduce corner and foot damage!
Easy to assemble corners for baskets for use during sterilisation, eliminating the risk of wrap perforation, caused by sharp corners or feet.

Additionally, the crepe version helps to reduce wet problems around basket corners.

Simply Fold, Tuck and Pinch to create a corner for baskets.

- Available in a heavy duty white rigid card.
- Fully Autoclavable.





Fold



Tuck



Pinch



CODE
CC001

DESCRIPTION
CLINIPAK WHITE RIGID CARD CORNERS

