

ERKA PRESSURE INFUSION

Pressure infusion devices are an essential item in any hospital. Pressure cuffs in combination with a manometer and inflation system provide for optimal pressure infusions. Suitable for solutions or blood bags in the following sizes: 500 ml and 1000 ml.



Features

- Sturdy nylon cuff in blue, with suspension loop for infusion stands
- Washable up to 60°C
- Precise reading of the level in the bag thanks to a transparent net on the cuff
- Manometer with 39 mm diameter aluminium scale up to 750 mm Hg

CODE

EP821 000 00
EP831 000 00

DESCRIPTION

ERKA PRESSURE INFUSER 500 ML
ERKA PRESSURE INFUSER 1000 ML

E&OE

18AUG08

