



ADULT CUFF & BLADDER SETS

Features

- Index and range markings
- Adult Range: 29-40cm
- Artery indicator label
- Latex Free inflation system
- Velcro fastener
- Extra Long Single Tube 80cm in length
- Available in Violet or Royal Blue



CODE	DESCRIPTION
ESCBVLF	LIBERTY CUFF & BLADDER SET VIOLET LATEX FREE 1 TUBE
ESCBARLF	LIBERTY CUFF & BLADDER SET ROYAL BLUE LATEX FREE 1 TUBE

