



## SPEY KITS DOG AND CAT

The Spey Kit comprises, amongst other important instruments, the HIFA Spey Hook. This is a long, thin, probe-like instrument with a curved end and a buttoned tip, used to hook around the uterus of a companion animal and to bring it into the incision where an ovariohysterectomy can be performed. Used as a means of neutering a female by removing the ovaries. The common 'spey' operation in dogs and cats.



SEE OVER FOR SPEY KIT CONTENTS AND ORDERING INFORMATION

E&OE First Published: 3 October 2008  
Last Modified: 8 January 2010  
File Name: info\_NKSC\_INSTRUMENTS\_SPEY-KITS





### SPEY KIT - CAT comprising:

- 2 Artery Forcep Mosquito curved 12.5cm
- 2 Artery Forcep Mosquito straight 12.5cm
- 1 Liberty Rochester Pean 14cm curved artery
- 1 Liberty Rochester Pean 14cm straight artery
- 1 Liberty Dressing Forcep 13cm blunt
- 1 Liberty Olsen Hegar 14cm needle holder
- 1 Liberty Mayo 15cm straight scissor
- 1 Liberty Surgical 13cm sharp/sharp straight scissor
- 1 Liberty Treeves Tissue Forcep 13cm 1x2T
- 1 Scalpel Handles No. 4
- 2 Bachus Towel clips 9cm
- 2 Allis Tissue Forcep 15.5cm 3/4T
- 1 Hifa Spey Hook with ball 14cm

Order Code: NKPROVSKC

### SPEY KIT - DOG LARGE comprising:

- 1 Artery Forcep Mosquito curved 12.5cm
- 1 Artery Forcep Mosquito straight 12.5cm
- 1 Liberty Rochester Pean 14cm curved artery
- 1 Liberty Rochester Pean 16cm curved artery
- 1 Liberty Rochester Pean 18cm curved artery
- 1 Liberty Rochester Pean 14cm straight artery
- 1 Liberty Rochester Pean 16cm straight artery
- 1 Liberty Rochester Pean 18cm straight artery
- 1 Liberty Dressing Forcep 16cm blunt
- 1 Liberty Olsen Hegar 14cm Needle Holder
- 1 Liberty Mayo 17cm straight scissor
- 1 Liberty Surgical 18cm sharp/sharp straight scissor
- 1 Liberty Treeves Tissue Forcep 16cm 1x2T
- 1 Scalpel Handles No. 4
- 4 Bachus Towel Clips 9cm
- 2 Allis Tissue Forcep 15.5cm 3/4T
- 1 Hifa Spey Hook with ball 21cm

Order Code: NKPROVSKDL

### SPEY KIT - DOG MEDIUM comprising:

- 1 Artery Forcep Mosquito curved 12.5cm
- 1 Artery Forcep Mosquito straight 12.5cm
- 1 Liberty Rochester Pean 14cm curved artery
- 1 Liberty Rochester Pean 16cm curved artery
- 1 Liberty Rochester Pean 14cm straight artery
- 1 Liberty Rochester Pean 16cm straight artery
- 1 Liberty Dressing Forcep 16cm blunt
- 1 Liberty Olsen Hegar 14cm Needle Holder
- 1 Liberty Mayo 15cm straight scissor
- 1 Liberty Surgical 16cm sharp/sharp straight scissor
- 1 Liberty Treeves Tissue Forcep 13cm 1x2T
- 1 Scalpel Handles No. 4
- 4 Bachus Towel Clips 9cm
- 2 Allis Tissue Forcep 15.5cm 3/4T
- 1 Hifa Spey Hook with ball 21cm

Order Code: NKPROVSKDM

### SPEY KIT - DOG SMALL comprising:

- 2 Artery Forcep Mosquito curved 12.5cm
- 2 Artery Forcep Mosquito straight 12.5cm
- 1 Liberty Rochester Pean 14cm curved artery
- 1 Liberty Rochester Pean 14cm straight artery
- 1 Liberty Dressing Forcep 13cm blunt
- 1 Liberty Olsen Hegar 14cm Needle Holder
- 1 Liberty Mayo 15cm straight scissor
- 1 Liberty Surgical 13cm sharp/sharp straight scissor
- 1 Liberty Treeves Tissue Forcep 13cm 1x2T
- 1 Scalpel Handles No. 4
- 4 Bachus Towel Clips 9cm
- 2 Allis Tissue Forcep 15.5cm 3/4T
- 1 Hifa Spey Hook with ball 14cm

Order Code: NKPROVSKDS

