



## HOT & COLD PACKS INSTANT COLD COMPRESS - UREA

The Liberty Instant Cold Compress temporarily relieves minor pain and swelling for sprains, aches and sore joints. It is conveniently disposable with no pre-chilling required for quick, effective relief. Product becomes cold on its own after activation.

The Liberty Instant Cold Compress is manufactured with urea instead of ammonium nitrate. Not only is this type of instant ice safer, but it is better for the environment. Convenient and portable, for use whenever and wherever cold therapy is needed. Store in a cool, dry place.

### Usage

Simply squeeze the pouch. The contents become immediately cold and remain cold for 14-40 minutes and achieves temperatures as low as  $-2^{\circ}\text{C}$

### Details INTICEUREA

Size: H228mm x W152mm  
Weight: 260g  
Pack Size: 40 per carton  
INTICEUREA INSTANT COLD COMPRESS (BOX 40)

### Details INTICEUREABOX

Size: H228mm x W152mm Individually boxed  
Box Size: H85mm x W115mm x D70mm  
Pack Size: 40 per carton  
INTICEUREABOX INSTAND COLD COMPRESS BOXED (CTN 40)

**Liberty**<sup>®</sup>  
HEALTH PRODUCTS



INTICEUREA

INTICEUREABOX

## REUSEABLE HOT/COLD GEL PACK

### Usage

To heat, microwave or place in hot water, follow instructions on the back of the compress. To obtain cold therapy, freeze for at least 1 hour.

This gel pack can reach a high of  $60^{\circ}\text{C}$  and a low of  $-12^{\circ}\text{C}$ . Maintains therapy for approximately 30 minutes.

### Details

Size: H125mm x W250mm  
Bag Style: Cloth feel  
Pack Size: 30 per carton  
INTGEL3510 REUSABLE HOT/COLD GEL PACK



INTGEL3510

E&OE First Published: 1 September 2009  
Last Modified: 8 October 2010  
File Name: info\_INT\_GEL-PACKS-HOT-COLD





## HOT & COLD PACKS REUSEABLE HOT/COLD GEL PACK

### Usage

To heat, microwave or place in hot water, follow instructions on the back of the compress. To obtain cold therapy, freeze for at least 1 hour. Can reach a high of 60°C and a low of -12°C. Maintains therapy for approximately 30 minutes. Can be used for either cold or hot. Remains soft and pliable after freezing. Cleans easily with soap and water.

### Details

Size:	H175mm x W125mm
Weight:	110g
Pack Size:	100 per carton
INTGEL3057	REUSABLE HOT/COLD GEL PACK



INTGEL3057

## BODY WARMER SELF ADHESIVE

### Usage

Open the outer package and shake gently, leaving the internal pack to activate with exposure to the air. It will gradually get warm and once temperature is reached it will last for around 12-15 hours, providing warmth for the user. Temperature can reach 50-68°C. Product provides long lasting warm to assist user in cold weather. Comes with adhesive to apply to clothing.

### Details

Size:	H175mm x W125mm
Pack Size:	240 per carton
INTBODYWARMER	BODY WARMER SELF ADHESIVE



INTBODYWARMER

## MAGIC CAT BODY WARMER

### Usage

To activate warmer, click metal disk by bending. Warmer has been activated if clear liquid begins turning white and solidifying resulting in heat. To reuse, place warmer in simmering water only for several minutes until all white crystals have dissolved. Remove from water and allow to cool. Your warmer is now ready to reuse.

### Details

Size:	H145mm x W110mm
Pack Size:	50 per carton
INTMAG4046CAT	INTCO MAGIC WARMER SELF ADHESIVE - CAT SHAPED



INTMAG4046CAT

