



MCCOY F.O. LARYNGOSCOPE FOR DIFFICULT INTUBATION

The McCoy Laryngoscope provides the clinician with greater flexibility and improved control in a wide selection of difficult intubation cases such as:

- Forward displacement of the Larynx
- Forward or prominent upper teeth
- Backward displacement of the tongue
- Decreased neck movement
- Cervical spine injuries
- Decreased mouth opening
- Recessive mandible

The McCoy Laryngoscope has a Macintosh style blade with a tip which is adjustable through 70°. A lever controls the tip angle during intubation to lift the epiglottis, giving a clear view of the vocal cords. Contact with the upper row of teeth and unnecessary pressure on the tongue and soft tissue are avoided.

Features:

- Swivelling blade tip facilitates difficult intubation
- Reduced likelihood of damage to teeth and soft tissue
- High quality stainless steel construction
- Autoclavable
- Fits all standard F.O. handles (Green System)

Liberty[®]
HEALTH PRODUCTS



Handle



McCoy Blade No. 3



McCoy Blade No. 4

CODE	DESCRIPTION
LLBMC3	MCCOY LARYNGOSCOPE BLADE NO. 3
LLBMC4	MCCOY LARYNGOSCOPE BLADE NO. 4
LLFH	LARYNGOSCOPE HANDLE
LLFB	LARYNGOSCOPE GLOBE

E&OE First Published: 19 February 2010
Last Modified: 19 February 2010
File Name: info_LL BMC_LARYNGOSCOPES_MCCOY

