



LARYNGOSCOPES FIBRE OPTIC

Features

- Autoclavable to 134°C
- Sure grip handle
- Takes 2 x 'C' size batteries (not included)
- Satin non-glare finish
- Cold fibre optic light
- Fully compatible with all fibre optic hook-on handles that conform to ASTM (F965)



CODE	DESCRIPTION
(NOTE: MCINTOSH BLADE IS CURVED)	
LLBFM0	MCINTOSH BLADE FIBRE OPTIC NO. 0
LLBFM1	MCINTOSH BLADE FIBRE OPTIC NO. 1
LLBFM2	MCINTOSH BLADE FIBRE OPTIC NO. 2
LLBFM3	MCINTOSH BLADE FIBRE OPTIC NO. 3
LLBFM4	MCINTOSH BLADE FIBRE OPTIC NO. 4
LLBFM5	MCINTOSH BLADE FIBRE OPTIC NO. 5
LLSFM3	LARYNGOSCOPE SET WITH 3 BLADES AND HANDLE IN CARRY CASE
(NOTE: MILLER BLADE IS STRAIGHT - NOT PICTURED)	
LLBFMI0	MILLER BLADE FIBRE OPTIC NO. 0
LLBFMI1	MILLER BLADE FIBRE OPTIC NO. 1
LLBFMI2	MILLER BLADE FIBRE OPTIC NO. 2
LLBFMI3	MILLER BLADE FIBRE OPTIC NO. 3
LLBFMI4	MILLER BLADE FIBRE OPTIC NO. 4
LLFH	FIBRE OPTIC HANDLE
LLFHAA	FIBRE OPTIC HANDLE FOR AA SIZE BATTERIES
LLHFS	FIBRE OPTIC HANDLE STUBBY
LLFB	HALOGEN BULB FOR FIBRE OPTIC HANDLE
LLB	REPLACEMENT OPTICAL BUNDLE
LLHFOA	LIGHT ASSEMBLY TO SUIT LLFH

E&OE First Published: 19 February 2010
Last Modified: 19 February 2010
File Name: info_LLFO_LARYNGOSCOPES_FIBRE-OPTIC

