




LARYNGOSCOPE FIBRE OPTIC SET

Features

- Stainless steel
- Autoclavable to 134°C
- Cold fibre optic light
- MacIntosh Blades in three (3) sizes: Size 1, Size 2 and Size 3
- Fully compatible with all fibre optic hook-on handles that conform to ASTM (F965)
- Halogen bulb
- Sure grip handle
- Satin non-glare finish
- Takes 2 x 'C' size batteries (not included)
- Carry case

Liberty[®]
HEALTH PRODUCTS



CODE
LLSFM3

DESCRIPTION
LARYNGOSCOPE SET FIBRE-OPTIC WITH 3 BLADES & CARRY CASE

E&OE First Published: 19 August 2008
Last Modified: 19 February 2010
File Name: info_LLSFM3_LARYNGOSCOPE-SET_FIBRE-OPTIC

