



## LARYNGOSCOPES STANDARD

### Features

- Autoclavable to 134°C
- Stainless steel handle
- Available in C or AAA batteries (not included)
- Satin Finish
- Fully compatible with all standard hook-on handles that conform to ASTM(F965)



CODE	DESCRIPTION
(NOTE: MCINTOSH BLADE IS CURVED)	
LLMS4	MCINTOSH SET 4 BLADES / HANDLE
LLBM0	MCINTOSH BLADE NO. 0 MATT
LLBM1	MCINTOSH BLADE NO. 1 MATT
LLBM2	MCINTOSH BLADE NO. 2 MATT
LLBM3	MCINTOSH BLADE NO. 3 MATT
LLBM4	MCINTOSH BLADE NO. 4 MATT
LLBM5	MCINTOSH BLADE NO. 5 MATT
BLBM1DISP	DISPOSABLE LARYN BLADE S/S MCINTOSH NO. 1
BLBM2DISP	DISPOSABLE LARYN BLADE S/S MCINTOSH NO. 2
BLBM3DISP	DISPOSABLE LARYN BLADE S/S MCINTOSH NO. 3
BLBM4DISP	DISPOSABLE LARYN BLADE S/S MCINTOSH NO. 4
(NOTE: MILLER BLADE IS STRAIGHT - NOT PICTURED)	
LLBMIM00	MILLER BLADE NO. 00 MATT
LLBMIM0	MILLER BLADE NO. 0 MATT
LLBMIM1	MILLER BLADE NO. 1 MATT
LLBMIM2	MILLER BLADE NO. 2 MATT
LLBMIM3	MILLER BLADE NO. 3 MATT
LLBMIM4	MILLER BLADE NO. 4 MATT
LLH	HANDLE FOR 2 C SIZE BATTERIES
LLHS	STUBBY HANDLE 2 AA BATTERIES (SIDE BY SIDE)
LLHAA	HANDLE FOR 2 AA SIZE BATTERIES
GLA	GLOBE - MCINTOSH / MILLER
GLFC	GLOBE - MILLER 00, 0 & 1

E&OE First Published: 19 February 2010  
 Last Modified: 19 February 2010  
 File Name: info\_LLS\_LARYNGOSCOPES\_STANDARD

