



LIBERTY OPHTHALMOSCOPE

The Liberty Ophthalmoscope fits into your pocket!

It's easy to use, lightweight and perfect for examination of the eye and it's background. Take your Liberty Ophthalmoscope with you everywhere you go throughout the day. Available in black only.

Features

- Halogen light for true tissue colour and consistent, long-lasting illumination
- Apertures include large and small spot, semi-circle and red-free filter
- Dust protection for maintenance free life
- Rubber brow rest prevents scratching of eye-glasses
- "AA" polycarbonate battery handle
- Pocket clip with integral on/off slide switch for ambidextrous use
- Supplied in a hard protective carry case

Technical Data

Correcting Lense Range:

+ 1 | 2 | 3 | 4 | 6 | 8 | 10 | 15 | 20 | 0 Diopter

- 1 | 2 | 3 | 4 | 6 | 8 | 10 | 15 | 20 Diopter


Aperture Wheel: Large spot, small spot, hemi-spot, red-free filter

Bulb: 2.5V 1.6W Halogen lamp

Batteries: 2 x "AA" batteries
or 2 x "AA" NiMH Rechargeable batteries
(Batteries Not Included)



Apertures



Liberty Pocket Instruments

*A complete range of quality diagnostic instruments for your pocket!
Top Liberty quality at a reasonable price!*

CODE	DESCRIPTION
LOPM	LIBERTY POCKET OPHTHALMOSCOPE

E&OE First Published: 19 August 2008
Last Modified: 27 July 2010
File Name: info_LOPM_OPHTHALMOSCOPE

