



## POCKET DIAGNOSTIC SET OTOSCOPE MINI/OPHTHALMASCOPE POCKET DIAGNOSTIC SET

The Liberty Pocket Otoscope Ophthalmoscope Diagnostic Set is super convenient, portable and value-packed with strong metal pocket clips. An ophthalmoscope giving 6 diagnostic beams, and a closed head otoscope, complete with 4 extra speculae, use both in the consulting room and when on the move. Supplied in a hard protective carry case



### Mini Fibre Optic Otoscope Features

- The Liberty Fibre Optic Mini Otoscope is a New and Improved Design, with a new switch and new modified bayonet fitting for a more secure attachment of speculae
- When the 2.5x Magnifying Glass is removed from the otoscope, external instruments can be safely inserted
- Otoscope has a bright light 2.5volt F.O. light source and battery handle (2 x AA size batteries)

### Ophthalmoscope Features

- Halogen light for true tissue colour and consistent, long-lasting illumination
- Apertures include large and small spot, semi-circle and red-free filter
- Rubber brow rest prevents scratching of eye-glasses
- Dust protection for a maintenance free life
- Pocket clip with integral on/off slide switch for ambidextrous use and battery handle (2 x AA size batteries)



| CODE      | DESCRIPTION   |
|-----------|---|
| LOTOPSET  | LIBERTY OTOSCOPE MINI/OPHTHALMASCOPE DIAGNOSTIC SET |
| LOTETD2.5 | DISPOSABLE EAR TIPS 2.5MM (PK 50)                   |
| LOTETD4   | DISPOSABLE EAR TIPS 4MM (PK 50)                     |

E&OE First Published: 1 October 2009  
Last Modified: 27 July 2010  
File Name: info\_LOTOPSET\_LIBERTY-DIAGNOSTIC-SET

