



## CLASSIC TUNABLE STETHOSCOPE

Dual head professional stethoscope with superior acoustic sensitivity and pressure sensitive diaphragm. The high performance stainless steel chest piece is at the top of its class. The definitive doctors stethoscope.



### New and Innovative Stethoscope

Light non-corrosive titanium design features "Superflex" binaural tubing assembly, colour coordinated plastic retaining, and non-chill rings. Complete with the innovative pressure sensitive diaphragm and advanced ear tip.

### Pressure Sensitive Diaphragm

This new design offers a pressure sensitive diagnosis, allowing clearer diagnosis of both high and low frequency sounds without having to change sides, as with standard stethoscopes.

### Advanced Soft Eartip

With its patented two-layered design, this model is extra flexible for added comfort.

### Binaural

"Superflex" deluxe stainless inner spring binaural assembly with ear tip.



### Features

- Tunable diaphragm
- Innovative design
- Optimum acoustic performance
- Finest craft finish
- Multi patents
- Soft ear tips
- Chest piece with non-chill ring
- Latex Free
- 5 year warranty
- Available in 5 colours
- Free name tag with each stethoscope
- Free spare ear tips and diaphragm
- Professionally packaged in a cardboard box





## CLASSIC TUNABLE STETHOSCOPE



### Colours Available

- Black**
- Raspberry**
- Navy Blue**
- Royal Blue**
- Purple**

- Colours are an indication only and may vary from actual product
- Colours may vary between supplies

### Technical Details

TYPE	Classic Tunable Dual Head Stethoscope	
MATERIAL	Head Piece	: Non-corrosive titanium
	Chest Piece	: Stainless steel brush finished
	Tubing	: PVC
	Ear Tips	: Latex free
	Diaphragm	: Hypo-allergenic insulating polymer
DIMENSIONS	Length	: 77mm
WEIGHT	Complete Unit	: 147gms
AUSTRALIAN DISTRIBUTOR	LIBERTY Health	
WARRANTY	5 Year	



CODE	DESCRIPTION
LSCLTB	LIBERTY CLASSIC STETHOSCOPE TUNABLE DIAPHRAGM BLACK
LSCLTNB	LIBERTY CLASSIC STETHOSCOPE TUNABLE DIAPHRAGM NAVY BLUE
LSCLTPU	LIBERTY CLASSIC STETHOSCOPE TUNABLE DIAPHRAGM PURPLE
LSCLTR	LIBERTY CLASSIC STETHOSCOPE TUNABLE DIAPHRAGM RASPBERRY
LSCLTRB	LIBERTY CLASSIC STETHOSCOPE TUNABLE DIAPHRAGM ROYAL BLUE

E&OE First Published: 14 November 2009  
 Last Modified: 25 February 2010  
 File Name: info\_LSCLT\_STETHOSCOPE-CLASSIC-TUNABLE

