



ELECTRIC PLASTER CUTTING SAW 125W INPUT POWER AND ALLOY BLADE

Liberty®
HEALTH PRODUCTS



Features

- Swing: Frequency: 133Hz
Range: ≥ 4 degrees
- Power Supply: 220V AC, 50Hz
- Input Power: 125W
- Blade Diameter: 62MM
- Power Cord Length: 2M
- Low noise total optimized mechanical structure and parts design < 85 dB
- Alloy blade is suitable for constant use
- Function of parenchyma protection is safe and reliable
- Variable speed control
- Suitable for cutting macromolecule plastic material



Complete With

- 1 Liberty® Electric Plaster Cutting Saw
- 2 Alloy Cutting Blades 62MM
- 1 Allen Key
- 1 Lockable Attache Style Carry Case

CODE	DESCRIPTION
PSEL	LIBERTY PLASTER SAW ELECTRIC
PSELB	LIBERTY PLASTER SAW SPARE BLADE

E&OE First Published: 4 September 2009
Last Modified: 9 February 2011
File Name: info_PSEL_PLASTER-SAW-ELECTRIC

