



7E-D PORTABLE SUCTION PUMP MULTI-VOLTAGE 15 LITRE



**The 7E-D is a
MULTI-VOLTAGE
unit and FULLY
PORTABLE**



Features

- Multi-Voltage
- Fully Portable
- Oil-free lubrication pump
- Square negative pressure dial
- Built-in bottle cavity
- Low noise
- Small size (H22CM x L29CM x D20CM)
- Easy to carry
- Fully portable unit incorporating three modes of electricity: 240V AC, 12V DC and battery
- Suitable for vehicles, clinic, home and outdoor places

Technical Information

Voltage:	AC 220V/240V 50HZ 12VDC
Power:	90 VA
Max Negative Pressure Value:	≥0.075 MPa
Pumping Rate:	≥15 Litres/Minute
Reservoir Volume:	1,000 mL
Noise:	≤65 Db (A)

CODE	DESCRIPTION
SPBB	LIBERTY 7E-D SUCTION PUMP 240V/12V/BATTERY 15L

SPARES/ACCESSORIES

SPBF	FILTER
SPBC12	SUCTION CANNULA 12FG
SPBSB1L	SPARE BOTTLE 1L

E&OE First Published: 14 October 2008
Last Modified: 14 May 2010
File Name: info_SPBB_SUCTION-PUMP

