



ANEROID SPHYG WITH TROLLEY & BASKET

Sphygmomanometer Features:

- Large, easy to read 150mm dial
- 300 mm Hg no pin stop mechanism
- Rear bracket storage compartment area
- Full 110° rotation for easy viewing
- Nylon velcro cuff, bladder, bulb and valve included
- Nickel plated air release valve
- 8 foot (when extended) of coil tubing

Trolley, Basket Features:

- Basket is a powder coated hard wire mesh
- Stainless steel pole
- Easy flow, high quality castors
- Broad base for extra stability
- Base options in light weight, medium weight and heavy weight



Powder coated hard wire mesh basket



Large, easy to read 150mm dial

Base Option Codes explained:
 LW = Plastic Base, Plastic Wheels (pictured)
 MW = Plastic Base, Metal Wheels
 HW = Metal Base, Metal Wheels



Broad base for extra stability

CODE	DESCRIPTION
STMALW	ANEROID SPHYG WITH TROLLEY & BASKET LIGHT WEIGHT BASE
STMALWF	ANEROID SPHYG WITH TROLLEY & BASKET LIGHT WEIGHT BASE LATEX FREE
STMAMW	ANEROID SPHYG WITH TROLLEY & BASKET MEDIUM WEIGHT BASE
STMAMWF	ANEROID SPHYG WITH TROLLEY & BASKET MEDIUM WEIGHT BASE LATEX FREE
STMAHW	ANEROID SPHYG WITH TROLLEY & BASKET HEAVY WEIGHT BASE
STMAHWLF	ANEROID SPHYG WITH TROLLEY & BASKET HEAVY WEIGHT BASE LATEX FREE
SCT	COILED TUBING FOR SPHYG (2 METRES) (NOT PICTURED)

E&OE First Published: 1 September 2009
 Last Modified: 28 April 2010
 File Name: info_STMA_ANEROID-SPHYG

