



POCKETSCOPE™ DIAGNOSTIC SET WELCH ALLYN® POCKET PROFESSIONAL

WelchAllyn®

Lightweight, compact diagnostic instruments offering high performance in a smaller package.

The Welch Allyn PocketScope™ uses halogen illumination. Built with many of the features found on 3.5V sets, it is smaller in size, but not in utility. The set includes otoscope and ophthalmoscope.

Features include:

- Two 2.5V handles (takes disposable AA batteries, included)
- Protective carry case
- Slick new ergonomic design
- Lens selection disc: offers choice of 48 lens settings
- Continuous aperture selection dial: offers six apertures (on patient side)
- Sealed optics: for improved long-term performance
- Choice of specula: with reusable or disposable Universal KleenSpec® specula
- Insufflator port: for pneumatic otoscopy
- Throat illuminator: built-in
- On/Off switch and rheostat control: let you set intensity of illumination



Bright, white light with Halogen HPX™

Fibre optics transmit light through a 360° bundle of optical fibres for accurate, shadow-free exams. The fibre optics extend to the tip of the instrument for more light output, less glare, and easier cleaning.

The Halogen HPX™ lamp system is an innovative lighting technology that uses a high-pressure mixture of halogen and xenon gases to provide 30% more light output than before.

What is a diagnostic set?

A basic diagnostic set consists of an ophthalmoscope for examination of the eye, an otoscope for examination of the ear, and a battery handle to power these instruments.

How do you use the throat illuminator on the PocketScope otoscope?

The built-in halogen throat illuminator offers the convenience of two instruments in one. Unscrew the otoscope head to reveal the high-intensity throat illuminator built into the otoscope base.

CODE	DESCRIPTION
WA728361	WELCH ALLYN POCKETSCOPE DIAGNOSTIC SET WITH CARRY CASE





POCKETSCOPE™ WELCH ALLYN® POCKET PROFESSIONAL

Improved Design! Improved Size and Shape! Superior Ergonomics!

Pocket Professional Otoscope

Key Features and Benefits

- Provide the advantage of halogen light and a fibre optic light transmission system
- Bright white cool light
- Ergonomically-designed handle with secure grip
- Attractive, modern appearance
- Sealed system with insufflator port for pneumatic otoscopy
- Accepts disposable KleenSpec (or reusable SofSpec®) specula
- Easy-to-use on/off switch
- Comes with 4 disposable specula
- Sliding back viewing lens



WelchAllyn®

Pocket Professional Ophthalmoscope

Key Features and Benefits

- Halogen Light coupled with precision optics cool, white, bright light
- Lens selection disc offers choice of 48 viewing-lens combinations
- Rubber brow rest prevents scratching of spectacles
- Sealed optics for improved long-term performance
- Continuous aperture selection dial offers 6 apertures (micro, small & large spot, fixation, slit and cobalt blue)
- Red filter across all apertures for easy identification of vessels and nerves fibres
- Cobalt blue filter can be used with fluorescein dye to quickly and easy detect corneal abrasions and foreign bodies
- Illuminated lens indicator makes lenses easier to select



CODE	DESCRIPTION
WA728353	WELCH ALLYN POCKETSCOPE OTOSCOPE WITH CARRY CASE
WA728357	WELCH ALLYN POCKETSCOPE OPHTHALMOSCOPE WITH CARRY CASE

